











The technical scheme do not consider the dimensions of any postural and anti decubitus cushions.

SW: Considered external tube to external tube of the frame. The dimension between the internal faces of the side guards is +10 mm than the SW.

SD: Considered the depth of the seat. It considers the dimension from the end of the seat to the back tube. SL: Considered the length of the seat pan. SL with CARBTS = SD - 20mm

CG: Considered from the backrest tube to the axle of the rear wheels.

SHF and SHR: Seat height front and rear can be adjusted only by tilting the seat. The frame tubes height to the floor is fixed.

BRH: Considered from the seat (upholstery or carbon) to the top of the backrest tube.

TOLERANCE: The product is manufactured with a tolerance of ± 5mm

For a right sizing of the wheelchair is recommended:

1 SW | Seat Width _ A space of ± 5 mm between the user basin and the side guards

2 SD | Seat Depth _ A space of ± 50 mm between the user's popliteal fossa and the rigid carbon seat.



SW: 340 - 460 mm in increments of 20 mm



SD: 380 - 480 mm in increments of 20 mm



BRH: 270 - 420 mm in increments of 15 mm



SHR: 360 - 430 mm in increments of 10 mm



SHF: 450 - 520 mm in increments of 10 mm



KHL: 300 - 440 mm in increments of 10 mm



CG: 100 - 150 mm in increments of 10 mm



BRA: 87° - 93° in increments of 3°

SDAL2 / BKRTT



CRW:



StF: 120 - 160 mm in increments of 20 mm



TOTAL WIDTH (TW) (Camber 0°) SW +175 mm (Camber 3°) SW + 220 mm



TOTAL DEPTH (TD) StF + SD + 105 + 115+ (100-CG) and (+ 50mm if with FRM93) and (- 20mm if with LWDF)



TRASPORT WEIGHT: (Without rear wheels) ± 6.0 kg*



MAXIMUM **USER WEIGHT:**



* CONFIGURATION: SW340 / SD440 / STF120 / LWDF / BRH255 / BKRUR / BRH255 / AXLMG0CV / F00TA2 / BRKASL / FRKDS2 / FRWHP4 /